




Happier Living Everyday



DIGITAL ECG HOLTER
TRAK-48

BPL TRAK-48

- **Compact Size, Light Weight**
68mm x 53mm x 16mm. 42 grams (excluding battery)
- **Fully Programmable, High Fidelity, Research Grade Performance**
Sample rate up to 1024 samples per second / No data compression
- **12-Lead/3-Channel Two Work Modes Switchable by Changing Patient Cable**
10 wire patient cable (I,II,III,aVR,aVL,aVF,V1-V6) / 7 wire patient cable (3 channel) or 5 wire patient cable (3 channel)
- **Lower Power Consumption**
Powered by on AAA size Alkaline / NIMH Battery \geq 48 hrs
- **Removable SD Card**
2GB records up to 48 hours of the 12 - lead ECG at 1024 SPS
- **128 x 64 Dot Matrix LCD Display**
Provides on demand ECG view / Confirms programming instructions / Visual confirmation of patient hook up
- **Pacemaker Detection**
- **Two Optional Modes of ECG Data Transfer**
Connected directly with PC through USB connector embedded in the recorder / Removable SD card

CERTIFIED
ISO
13485
COMPANY