



**Transport ventilation
without compromise**

DRÄGER
OXYLOG® 3000 PLUS

Rising up to the challenge



ICU ventilation standard at all times

Whether you're in the field or in the hospital, successful treatment is the result of teamwork. Medical professionals work together to save lives as the patient progresses through a chain of events – from rescue and stabilization to diagnostics, from treatment to recovery and, finally, to discharge. However, a chain is only as strong as its weakest link, and for a variety of reasons, patients are most vulnerable during transport. Maintaining the highest possible standard of ventilation at all times, even during transport, can mean the difference between success and setback.



OXYLOG 3000 PLUS – READY, COMPACT, ROBUST

At Dräger, we took on the challenge – raising the bar in transport ventilation. Drawing on over a century of experience in the field of mechanical ventilation, we created a device that gives you state-of-the-art ventilation in a compact and rugged format that is simple and intuitive to operate. Now, you can confidently provide high-quality ventilation therapy for pediatric and adult patients during transport – anytime, anywhere.

OXYLOG® 3000 PLUS

Ventilation modes	VC-CMV, VC-AC, VC-SIMV, SpnCPAP, PC-BIPAP (PC-SIMV+)
Additional settings	NIV and AutoFlow, Apnea ventilation, Pressure support
Options*	AutoFlow, capnography, real time data export
Operating time (new battery, fully charged, typical ventilation, without CO ₂ sensor, reduced display brightness)	Appr. 9.5 hours
Operating time (new battery, fully charged, typical ventilation)	Appr. 7.5 hours
Weight	Appr. 5.8 kg (basic unit with internal battery)

* Options can be purchased during the initial ordering process or as future upgrades.

ICU ventilation performance for transport

HIGH PERFORMANCE VENTILATION CAPABILITY

- Wide range of pressure and volume controlled ventilation modes
- Adjustable inspiration termination and apnea backup ventilation
- Low oxygen consumption because no basic flow is needed

NON-INVASIVE VENTILATION

- Synchronization with spontaneously breathing patients
- Patient ventilation with automatic leakage compensation

AUTOFLOW® (OPTIONAL)

- ICU ventilation performance to minimize peak airway pressures
- Applicable with all volume controlled modes

INTEGRATED CAPNOGRAPHY (OPTIONAL)

- Mainstream CO₂ sensor for short response times with plug & play functionality
- Integrated monitoring: etCO₂ and waveforms

READY FOR PEDIATRIC PATIENTS

- Accurate tidal volume delivery to 50ml
- Dedicated pediatric breathing circuit



MT-5800-2008



ENHANCED MONITORING

- Uninterrupted displaying of numeric and curve data
- BTPS and automatic altitude compensation

DATA MANAGEMENT (OPTIONAL)

- Real-time data transfer
- Remote service to transfer Oxylog 3000 plus status to DrägerService

SIMPLE OPERATION

- Quick start time, no need of Flow Sensor calibration
- High-visibility display and intuitive user interface

PERFECTLY PORTABLE & PROVEN DURABILITY

- Simple but effective mounting systems for safe, hands-free operation
- Powerful battery with up to 9,5 hours of independent operation

WIDE RANGE OF ACCESSORIES

- Availability of standardized Dräger accessories for the Oxylog family
 - Wall holder, carrying solutions, trolleys and many more

Setting standards for transport ventilation



With the Oxylog 3000 plus, you can choose the system you need today while maintaining a variety of options for the future. A full range of pressure and volume controlled ventilation modes combined with NIV, pressure support, incremental O₂ concentration regulation and apnea ventilation enable high-quality therapy and provide confidence during transport.

With its advanced features and functions, the Oxylog 3000 plus sets standards for transport ventilation that will help you rise up to the challenges of modern transport ventilation – now more than ever.

- Intra- and inter hospital transport
- Comprehensive treatment for emergency care and transport
- Uninterrupted ICU ventilation therapy
- No compromise when transporting ventilated patients
- Enhanced and effective hospital workflow
- Two years service interval without the need of a display change

Accessories

D-18571-2010



NIV FULL-FACE MASK CLASSICSTAR®

Dräger offers a wide variety of dedicated accessories for the Oxylog 3000 plus. The ClassicStar NIV full-face mask is designed for non-invasive ventilation therapy. With its air-filled cushion, the mask can be adjusted individually to the patient's face and so provides effective therapy.

Education

9067073



DRÄGER ACADEMY – BASICS OF RESPIRATION AND VENTILATION

How exactly does breathing work? What happens when a person actually "runs out of air" and must be ventilated? All these and various other respiration and ventilation questions are answered by the Dräger eLearning program.

Training

9066314



OXYLOG 3000® PLUS PRODUCT TRAINER

Dräger supports your everyday clinical work with online training offers. With our interactive product trainer you get useful information about the Oxylog 3000 plus and you can easily practice the intuitive control concept and learn operating the device.

Service

D-32436-2011



DRÄGERSERVICE – BECAUSE QUALITY COUNTS

For any requirement you have, we have the service contract that fits: from inspection-only and preventive maintenance to complete service packages. We have the answer for your individual needs.

For further information and downloads please visit www.draeger.com

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

Locate your Regional Sales
Representative at:
www.draeger.com/contact



**REGION EUROPE CENTRAL
AND EUROPE NORTH**
Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH
Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA
Drägerwerk AG & Co. KGaA
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC
Dräger Médical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA
Dräger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free+1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

**REGION CENTRAL
AND SOUTH AMERICA**
Dräger Panama Comercial
S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com